## ENGLISH

DUO CONTACT LENS CASE

SMART COLOUR CAP SYSTEM

3D CHARACTER CONTACT LENS CASE

NO LEAK LOCKING SYSTEM

DUO CONTACT LENS CASE

- For Soft, Hard or Gas permeable lenses
- FDA and CE approved

### DIRECTIONS FOR USE

- Before initial use rinse lens case and caps thoroughly with fresh contact lens disinfecting solution. Use this method to clean the lens case weekly. Do not use soaps, detergents or other household cleaners.
- After you remove your lenses from the lens case, empty & rinse it with fresh contact lens disinfecting solution. Allow the case to dry thoroughly.
- Always store lenses in fresh solution. Never reuse solution.
- Follow your manufacturer's instructions for your lens care regime.

### Warnings:

- Lens cases can be a significant source of microbial contamination. To help prevent eye infections, lens cases should be properly cleaned, rinsed and air dried.
- Lens cases should be replaced every month, or as directed by your eye care practitioner.
- Use of this lens case with heat may cause warpage. **USE FOR STORAGE DURING CHEMICAL DISINFECTION ONLY. DO NOT USE WITH HEAT.**
- Some chemical & hydrogen peroxide systems require special lens cases. Follow the manufacturer's instructions or consult your eye care practitioners.

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